

Bio-Speedy vNAT Viral Nucleic Acid Buffer



Bio-Speedy®vNAT® Viral Nucleic Acid Buffer provides an innovative, rapid and effortless solution.

Bio-Speedy®vNAT® Viral Nucleic Acid Buffer enables viral nucleic acid extraction from nasopharyngeal aspirate and lavage, bronchoalveolar lavage, saliva, nasopharyngeal swab, oropharyngeal swab, and sputum samples.

After the samples are combined with Bio-Speedy®vNAT® Viral Nucleic Acid Buffer, within 5 minutes the vNAT® reagent breaks down the virus and releases the nucleic acids. Thus, kit allows going from sample to qPCR within 5 minutes.

Quantity: 100 ml

Box Size: 10cm x 8cm x 6.5cm

Catalog Number: BS-NA-510

